



A New Quantitative Method for Analysis of Drill Cuttings and Core for Geologic, Diagenetic and Reservoir Evaluation

Jon Sliwinski

SGS Minerals Services, Lakefield, ON, Canada

jonathon.sliwinski@sgs.com

and

Michael Le Strat and Murray Dublonko

Breaker Energy Ltd., Calgary, AB, Canada

Summary

Drill cuttings have traditionally been utilized as a first look, qualitative tool but are often not re-examined once wireline logs are run. The advent of advanced technology from the mining industry has broad applications for quantitative rock analysis of oil and gas reservoirs. The automated, high resolution data provided from QEMSCAN[®] analysis of drill cuttings provides much more information than conventional optical analysis.

SGS and Breaker Energy have completed a pilot study to evaluate the applications of this new technique on a recently drilled well. The horizontal well in NE British Columbia targeted the Doig Formation, a tight, shaly gas sand with localized phosphatic zones. The QEMSCAN[®] data provided mineralogical and textural information on the distribution of coarse grains and shale, visual confirmation of cementation mineralogy and fabric as well as quantitative porosity numbers. The analysis aided in developing the depositional model, diagenetic history, reservoir characteristics and completion/production strategies. This paper will present the methodology, data and conclusions of this study.

Introduction

Drill cuttings are typically the first piece of data available from the subsurface during drilling. These samples of the subsurface formations are collected at surface and used by wellsite geologists to prepare a striplog and then typically disappear from the data stream. These samples can provide valuable information and in the absence of core provide the only actual contact with the rock. Recently there has been renewed interest in the data available from drill cuttings and the application of advanced techniques (Scanning Electron Microscopy, X-ray Dispersive Energy Spectroscopy, Back Scattered Electron Imaging and Digital Imaging and Processing) as well as advanced sample preparation (grain mount thin sections).

QEMSCAN[®] is an automated mineralogy analysis system that has a long history in the mining sector. It is a combination of a Scanning Electron Microscope (SEM) and Electron Dispersive Spectroscopy (EDS) with automated analysis software that uses BSE (Backscattered Electron Intensity) and EDX (Energy Dispersive X-ray) spectra to identify the mineral or phase present at the analysis point. QEMSCAN[®] can measure and identify up to 600 points per second. This speed allows for detailed mineral mapping of samples, with the production of digital images, and superb statistical data on the modal mineralogy; typically millions of data points are measured during each analysis. QEMSCAN[®] is able to analyze cuttings, as well as conventional and rotary cores. Prepared samples are polished sections from cuttings, polished thin sections from cuttings or core, and solid core.

Method

Drill cuttings were provided from the wellsite and first examined and photographed using a traditional stereoscopic microscope (Figure 1). The wellsite geologist also provided his observations in a traditional striplog for comparison. Polished sections were prepared using micro-riffling samples from each five meter depth interval. The riffling process ensures that a random selection of particles from the interval is present in the analysis. These particles are then prepared in a resin polished section for analysis by QEMSCAN[®].




Figure 1: Stereo Microscope photographs of the cuttings and one sandstone particle from a common depth interval.

When analyzing samples using QEMSCAN[®], there are several choices for operation mode (Figure 3), resolution and magnification. There are several scanning modes depending on what information is sought, for the purposes of this study we ran a line scan (**Bulk Modal Analysis**) for overall mineralogy and a particle scan (**Particle Map Analysis**) based on the use of cuttings as opposed to core or thin section. In the mining industry, samples are usually crushed and screened into size fractions and the resolution and magnification is dependent on the size of the fraction being analyzed. In this case we did not screen the sample as the particle size was fairly consistent at 1-0.5 mm, which also determined the magnification of 20x in order to capture the particles in the view window. The matter of resolution was not as easily determined and several different resolutions were compared. In the end it was determined that a resolution of 5µm was optimal for showing details of the grains and giving a clear view of the porosity and cementing features. In order to save time in the analysis all particles less than 250 µm were not analyzed. As a control the modal mineralogy from the BMA, run on the entire sample, was compared to the mineralogy from the PMA to ensure the analysis was statistically accurate.




Figure 3: QEMSCAN modes of operation.

Data

Presented in Figure 4 is the data from one polished section prepared from a depth interval in the Doig Formation. The first image shows the total particles scanned. In this case there were 197 particles larger than the lower limit of 250µm.

The second image is one of the sandstone grains enlarged showing the detail. In the false colour map pink represents quartz, orange potassium feldspar, red pyrite and turquoise dolomite. The grain is very clearly cemented with a dolomitic cement as well the pore space, the white voids, are lined with dolomite suggesting the formation may respond well to an acid treatment prior to the planned frac. Porosity of this particle can be determined by the analysis software and combined with that of the other grains giving a quantitative value for the porosity.

The third image is the individual quartz and potassium feldspar grains from the sandstone which were then analyzed by the software for roundness (shape factor), and average grainsize. This textural information can then be incorporated in a depositional model.

The fourth image is the pore throats of the sandstone grain. The analysis software determined the porosity to be 5.06% and have an average size of 10.33µm. A new porosity value can then be calculated using a minimum “effective” pore size cutoff, giving a value for effective porosity. Alternately part of the dolomitic pore lining could be added to the porosity, simulating an acid stimulation, giving an “improved” porosity value.




Figure 4: a) False colour maps of the 197 particles mapped in the polished section made from the drill cuttings of the same depth interval as the photographs in Figure 1; b) false colour and BSE maps of one of the sandstone particles from the 197 cuttings c) The quartz and feldspar grains within that sandstone particle, and in black below the pore throats.

Conclusions

The QEMSCAN[®] provided useful data in the evaluation of this reservoir. The most unexpected feature was the dolomitic cement which provided valuable input into the completion strategy for this well. Additional data is still being incorporated by the geologists. The speed and capacity to evaluate large volumes of samples from this horizontal well in a small amount of time, allowed a detailed analysis of small changes within the reservoir.

QEMSCAN[®] demonstrated it is a valuable new tool to the oil and gas industry. It has broad applications from simple analysis of cuttings for modal mineralogy to detailed mapping of core. It has applications in the analysis of individual wells as well as showing potential for large field studies and field to field comparisons.

Acknowledgements

The authors would like to acknowledge SGS and Breaker Energy for providing the resources to complete this study.